

# INNER LOK™ ANKLE BRACE

Ideal for Sports to Support Weak or Injured Ankles.

**40-511 Series**

COLOR: BLACK

## The Ankle Brace with Inner Strength

### JUNIOR BALLISTIC NYLON CONSTRUCTION

Internal foam padding is laminated to a junior ballistic nylon exterior. The foam padding interior provides a comfortable fit while the junior ballistic nylon exterior makes the brace extremely durable. This combination of materials is also lightweight and fits easily in your regular athletic shoe.

### INTERNAL FIGURE EIGHT STRAPS

This revolutionary design has two figure eight straps that are actually pre-configured through the inside of the brace in the correct position to assure proper strap placement every time. The figure eight straps provide extra support and allow for a more custom fit. The straps are made of a heavy duty DuPont Cordura® material and are protected inside a layer of junior ballistic nylon for maximum durability.

### SHORT CURVED SEAMLESS ARCH

The short curved arch fits the contour of the foot and allows full range of motion comfortably. The seamless arch virtually eliminates irritation to the bottom of the foot.

### BREATHABLE TONGUE

This advanced tongue material wicks moisture away from the skin. At the same time, the Breathable Tongue allows air to pass through to the skin, keeping your ankles cooler and drier. The breathability of this fabric also helps eliminate skin irritation and moisture problems associated with other fabrics. The tongue is completely sewn in to stop migrating and prevent pinch points.

### FULL ELASTIC BACK

Only Swede-O has elastic all the way up the back of the brace to ensure complete unrestricted blood flow to the Achilles' tendon. The full elastic back virtually eliminates the chance for blistering.

### EXCLUSIVE CLOSE SPACED EYELETS

Three close spaced eyelets are located in the middle of the brace, directly opposite the heel opening, to lock the heel in place. This keeps the brace in the correct position and prevents slippage on your ankle.



### TO DETERMINE SIZE:

Size according to shoe size:

Size	Men's Shoe	Women's Shoe
Small	6 - 7	7 - 9
Medium	8 - 10	10 - 11
Large	11 - 12	12 - 13
X-Large	13 - 14	14+

### WEARING INSTRUCTIONS:

Wear brace over sock. Slip foot into brace. Tighten the InnerLOK8™ straps for a secure fit. First, adjust the inside strap until snug and repeat with the outside strap. Tighten laces from bottom to top. Retighten InnerLOK8™ straps for an even and comfortable fit.



Swede-O and InnerLOK 8 are Trademarks of Swede-O, Inc.

**800-327-4110**

**F·L·A ORTHOPEDICS®**

Website: [www.flaorthopedics.com](http://www.flaorthopedics.com)

©2002 FLA ORTHOPEDICS, INC.